SEFFICARE®

Autologous Tissue Treatment

EASY - SAFE - STANDARDIZED

USER MANUAL

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DEVICE DESCRIPTION AND USE

The SEFFICARE® kit is designed for **physicians** who want to tap in the potential of microfragmented adipose tissue in the self-repair of integumentary and musculoskeletal tissues.

Several studies have shown that the subcutaneous adipose tissue is made of adipose cells and stromal tissue (SVF) containing mesenchymal cells and especially adipose derived stem cells (ADSCs). Thank to the natural presence of these cells, the microfragmented adipose tissue enhances the self-repair potential of tissues.

Under this procedure the tissue harvested and prepared with SEFFICARE® is implanted in the same patient (autologous graft), in order to obtain the regeneration, reconstruction or self-repair of integumentary and musculoskeletal tissues.

This device can be used on different body areas on the same patient and within the same procedure.

SEFFICARE® is a disposable sterile kit to be used on a single patient.

SEFFICARE® is designed to facilitate and standardize the harvesting, preparation and grafting of subcutaneous adipose tissue.

ABSOLUTE CONTRAINDICATIONS

- Ongoing infections in the area of tissue harvesting or grafting.
- Presence of malignancies in the area of tissue harvesting or grafting.
- Pregnancy or breast-feeding.
- Anticoagulant therapy or severe coagulation disorder.
- Allergy to local anesthetic.
- Dysmorphophobia.
- Ongoing immunosuppressive therapies.
- Debilitated subjects.

RELATIVE CONTRAINDICATIONS

- Strong tobacco addiction.
- Presence of non-resorbable material in harvesting or grafting area.

WARNINGS

- 1. This is a disposable device.** Reusing the SEFFICARE® system may cause infections and the transmission of communicable diseases.
- 2. This device is sterile.** If the sterile packaging appears damaged, do not use the product.
- 3. This device cannot be re-sterilized.** Do not autoclave.
- 4. Do not use the device after the expiry date on the label.**
- Before use and during transport of the device, the **surrounding temperature must range between 5°C and 40°C.** Excessive temperatures may compromise the reliability of the device.
- Utmost caution is necessary when using this device on patients in precarious state of health.
- This device is not meant for removal of adipose tissue in a body remodeling treatment or weight reduction.
- During the procedure do not force the guide rod. In case of breakage of the rod or of other parts of the guide **interrupt immediately** the procedure.
- The harvesting procedure must be performed in body areas where the pinch test is = or > 4 cm.
- In order to guarantee the traceability of the product, please put the sticker identified by “**Consenso - Consent**” and **QR code** located on the bottom of SEFFICARE® box on the form for the informed consent to be signed by the patient.

CAUTIONS

- The use of this device is limited only to those physicians who have been specially trained.
- The procedure outcomes vary depending on the patient's age, state of health, life habits and surgical site as well as on the physician's experience.
- Results achieved may not be definitive.
- The quantity of tissue harvested must be limited to what is actually needed.
- The surgical field needs to be completely sterile.** Please follow medical standards designed to guarantee the maintenance of sterility during the procedure (e.g. do not touch the non-sterile parts of the patient or operator such as garments, jewellery, glasses, phones, etc.).

UNDESIRE EFFECTS

Undesired effects deriving from the procedure may be: asymmetries, infection, inflammation, hematomas, ecchymoses, skin irregularities, scar healing defects, skin dyschromia, fat necrosis, oil cysts, allergic reactions.

SEFFICARE
LIPOUSION SYSTEM
GMDN 56627

CE 0425

REF SEFFICARE

LOT 2004L3

05-05-2020

STERILE EO

2023-04

Consenso - Consent

Seffiline S.r.l.
Via Delle Lame, 98
40122 Bologna - Italy
info@seffiline.com
www.seffiline.com

(01)18053300130013(17)230430(10)2004L2

KIT CONTENTS SEFFICARE®

THE STERILE BAG INSIDE THE BOX CONTAINS:

- no. 1 guide
- no. 1 1.0mm-hole infiltration and harvesting cannula
- no. 1 multiple syringe decanting stand
- no. 1 siringa VacLok® 10 ml
- no. 3 10 ml luer lock syringes
- no. 3 3 ml luer lock syringes
- no. 3 plugs for luer lock syringes
- no. 2 connectors for luer lock syringes
- no. 1 cylindrical container
- no. 4 needles with different caliber and length
- no. 4 gauze swabs
- no. 1 drape
- no. 1 drape with opening

KIT UNPACKAGING

STEPS IN WHICH YOU DO NOT NEED TO WEAR GLOVES

(not present in the package).

1. Take the package out of the box keeping the label upwards.

2. Open the package by pulling apart the upper and lower sides of the package.

3. Slide the package content on the Mayo table, keeping the label upwards. Do not touch the inner content.

STEPS IN WHICH YOU MUST WEAR STERILE GLOVES

(not present in the package).

4. Wear sterile gloves.

5. Remove gently label and tape.

6. Grab the drape edge wrapping the kit and spread it on the Mayo table without touching the table top.

7. Repeat procedure no. 6 for the other 3 edges.

8. Take the drape with opening on the top of the tray and place it to the side of the tray.

9. Place the kit components on the drape as is more convenient for the operator.

STEPS IN THE PROCEDURE

The procedure is to be performed in sterile conditions: the skin should be carefully cleaned with a disinfecting solution both in the harvesting area and in the implantation area prior to the procedure. Do not compromise the sterility conditions during the procedure. You should absolutely avoid any contact with non-sterile parts of the patient or operator (e.g. do not touch the non-sterile parts of the patient or operator such as garments, jewellery, glasses, phones, etc.)

HARVESTING

The area where the tissue is harvested should be chosen based on the following criteria:

- adequate quantity of subcutaneous adipose tissue;
- size of the area: the harvesting area should be at least 20 cm in diameter;
- the most common harvesting areas are: abdomen, hips and the trochanteric region.

Harvesting phases:

A. PREPARATION

- A1.** Clean the skin with a disinfecting solution.
A2. Lay the drape with opening in position.

B. INFILTRATION OF LOCAL ANAESTHESIA

The local anaesthetic is NOT included with the SEFFICARE® kit. Using local anaesthetics containing adrenalyn is recommended.

B1. Choose the spot for insertion of the cannula: it should be 1 cm out of the perimeter of the harvesting area (about 20 cm in diameter).

B2. Use the 27G needle to inject the anaesthetic solution.

B3. Make an opening in the skin with the 18G needle.

B4. Connect the multi-hole harvesting cannula with the 10 ml luer lock syringe containing local anaesthetic.

B5. Insert the syringe with the cannula into the guide.

B6. The cannula must be inserted into the hilt of the guide and connected to the 10 ml VacLok®.

B7. Rotate the syringe so that the barrel flanges get caught in the locks at the back of the guide.

B8. Insert the tip of the cannula into the skin opening perpendicularly until the guide stops.

B9. Rotate the cannula by 90 degrees so as to set the guide parallel and adhering to the skin plane.

B10. Infiltrate the local anaesthesia with a back and forth movement of the cannula while maintaining the guide parallel and adhering to the skin. The anaesthetic solution should be injected in the whole harvesting area. If more anaesthetic is necessary, repeat the same procedure from step **B5** after recharging your syringe with more anaesthetic solution.

B11. WAIT 10 MINUTES.

C. HARVESTING OF TISSUE

C1. Connect the multi-hole harvesting cannula with the 10 ml VacLok® syringe. Make sure the VacLok® syringe graduated scale is turned upwards and visible.

C2. Insert the syringe with the cannula into the guide.

C3. The cannula must be inserted in the guide hilt.

C4. Rotate the syringe so that the barrel flanges get caught in the locks at the back of the guide.

C5. Insert the tip of the cannula into the skin opening perpendicularly until the guide stops.

C6. Rotate the cannula by 90 degrees so as to set the guide parallel to the skin plane and the cannula in the subcutaneous plane.

C7. While holding the guide with your dominant hand, pull the syringe plunger back and rotate it clockwise so as to block the plunger and create the vacuum in the syringe barrel.

C8. While maintaining the guide parallel and adhering to the skin, move the cannula back and forth with a fan movement in the whole harvesting area. The cannula must not come out of the skin hole beyond the reference line on the guide. The tissue will be sucked into the barrel. Stop as soon as you have collected 5 ml of tissue in the syringe.

C9. Pull out the cannula from the skin.

C10. Rotate the syringe until the plunger flange is released from the lock at the back of the guide.

C11. Pull out the syringe, disconnect the cannula and attach the transfer.

C12. Connet the 10ml syringe to the other opening of the transfer, then rotate the VacLok® plunger counter clockwise for ¼ of a turn and transfer the tissue into the 10 ml syringe.


C13. Fill the syringe with sterile physiological solution (NOT included in the kit) you previously poured into the cylindrical container.

C14. Use the luer lock plug to close the syringe.

C15. Put the syringe (plug down) in the multiple syringe decanting stand. Repeat the procedure from step C1 to step C15 until you have from one to three 10 ml syringes of tissue according to the treatment plan.

PREPARATION OF TISSUE

After some minutes the syringes in the decanting stand will show a separation by gravity of the tissue (top) from the physiological solution (bottom).

D. WASHING THE TISSUE

D1. Pull out the syringe from the multiple stand, hold it in a vertical position (DO NOT TILT nor TURN UPSIDE DOWN), remove the plug, discard the washing physiological solution while retaining all the tissue inside the syringe.


D2. Put the plug back in position and the syringe back in the stand.

D3. Repeat the steps **D1** and **D2** for all the syringes in the stand.

E. FLUIDIFICATION

E1. Pull out the 10 ml syringe containing the tissue from the stand and remove the plug.

E2. Use the connector to connect the syringe to an empty 10 ml syringe.

E3. After pushing the tissue about 10 times from one syringe to the other, a good fluidification of the tissue is obtained.

E4. Disconnect the syringes and use the plug to close the syringe with the tissue.

E5. Put the syringe in the stand.

E6. Repeat the procedure from step **E1** to step **E5** for all the syringes.

GRAFTING

F1. Remove the plug from the 3 ml syringe containing 2.5 ml of tissue.

F2. Disinfect thoroughly the grafting area.

F3. Proceed by injecting the tissue in the anatomical area and in the quantity requested by the treatment indication with a suitable needle.

DISPOSAL

In case of use on patients with contagious bacterial and/or viral infections, some components could result highly contaminated.

Once used, do not dispose of the device with domestic waste. The device must be treated as a **Special Waste**.

Please deliver your used device to an authorised collection center.

Contact the Local Authorities to learn about the disposal regulations of such waste.

In order to respect the environment and not to cause any damage to health, we ask

you to comply with these recommendations.

The presence of this symbol highlights the obligation of differentiated waste disposal.

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The product SEFFICARE® complies with the requirements of the directive 93/42/CEE on medical devices.

SEFFILINE

**Via delle Lame, 98
40122 Bologna - Italy
+39 051 0217454**